

## Examline 740 Gynecological Examination Table

Ref: M2S 300 0 0



Bıçakçılar Examline 740 Gynecological Examination Table is designed for gynecological, rectoscopy/proctology, urology and urodynamic examinations, minor surgical procedures and obstetric deliveries. With a soft-touch ergonomic design. Examline 740 is constructed with high quality materials providing effective and long operating life in hospitals and clinics.

Gynecological Examination Table consists of backrest and seat sections. Seat extension pads are offered as optional accessory. As standard accessory either foot supports or knee crutches integrated with handles can be selected according to customers' need.

All surfaces contacting with the patient are covered with hygienic polyurethane pads. Antistatic pads are waterproof and easy to clean.

Below listed tabletop positions of Examline 740 can be adjusted with remote control or with multi-functional foot pedals:

- Tabletop height adjustment,
- Trendelenburg/reverse Trendelenburg positions,
- Backrest section forward and backward movements,
- Handle integrated foot supports and knee crutches up and down movements.

Three positions can be preset in memory by use of the remote control unit. All positions are adjusted by electrically operated pistons. Gynecological Examination Table is equipped with battery recharge system. Movements of the table are fed from the battery circuit. Battery power has 30 minutes operating time (approx 30 full positioning).

Electrically adjustable foot supports and a low starting position (600 mm mounting height with 5° reverse Trendelenburg position) offer maximum comfort. The patient can take her place easily - without the need of foot step.

An examination tissue paper dispenser is stored under the backrest section.

A stainless steel waste tray is located below the seat section.

Examline 740 is available in blue color.

## Technical Details

### Height adjustment:

660 - 910 mm

### Tabletop: 550 x 1270 mm

### Table base:

600 x 1085 mm

### Head section:

550 x 320 mm

### Backrest section:

550 x 600 mm

### Seat section:

550 x 385 mm

### Length of seat section

w/seat extension pad:

550 x 550 mm

### Voltage:

230V, 50/60Hz

## Packing Details

### Net weight (w/o accessories):

180 kg

### Gross weight (w/o accessories):

225 kg

### Packing dimensions:

1000 x 1680 x 1100 mm

## Positions

### Backrest section:

52° upward from horizontal position

### Seat section:

In combination with backrest section, 19° Trendelenburg and 5° reverse Trendelenburg (providing 600 mm mounting height)

### Foot supports:

Integrated with handles, can be swiveled to the sides, and moved up and down (to 140 mm from the floor)

### Knee crutches:

Integrated with handles, can be swiveled to the sides, and moved up and down.



## Accessories For Examline 740 Gynecological Examination Table

### Standard Accessories

- ▶ Knee crutches (2 ea)
- ▶ Waste tray ( 1 ea)
- ▶ Examination tissue paper (1 ea)



**Knee crutches**  
Ref: M2S 325 5 0 (Ref of 1 ea)



**Waste tray**  
Ref: M2S 171 1 0



**Examination tissue paper**  
Ref: MT1 06 0001

### Accessories

- ▶ Foot support
- ▶ Anaesthetist's frame
- ▶ I.V. stand
- ▶ Arm boards
- ▶ Seat extension pad
- ▶ Leg section
- ▶ Examination light (Luxline 40 LED)
- ▶ Physician stool
- ▶ Physician stool w/foot control
- ▶ Physician chair
- ▶ Physician chair w/foot control



**Foot support**  
Ref: M2S 335 5 0 (Ref of 1 pair)



**Anaesthetist's frame**  
Ref: M3S 000 0 0



**I.V. stand**  
Ref: M3T 000 0 0



**Armboards**  
Ref: M9F 200 0 0 (Ref of 1 ea)



**Seat extension pads**  
Ref: M2S 340 0 0 (Ref of 1 pair)



**Leg sections**  
Ref: M2S 360 0 0 (Ref of 1 pair)



**Examination light (Luxline 40 LED)**  
Ref: M1V 400 0 0



**Physician stool**  
Ref: MOB 000 0 0



**Physician stool w/foot control**  
Ref: MOB 200 0 0



**Physician chair**  
Ref: MOB 100 0 0



**Physician chair w/foot control**  
Ref: MOB 300 0 0