

## Surgiline 2000E Surgical Table

Electrohydraulic

Ref: M2C 650 0 0



Biçaklılar Surgiline 2000E Surgical Table is designed to accommodate different patient positions required for demanding surgical procedures, with an ultimate focus on patient and user safety. Surgiline 2000E Surgical Table is constructed with stainless steel and high quality materials to provide durability and long operating life.

Surgical positions of Surgiline 2000E are adjusted by an electrohydraulic system. Surgiline 2000E Surgical Table has up/down, sliding, Trendelenburg/reverse Trendelenburg, right/left lateral, renal surgery, Trendelenburg and back section upward combination, and back section tabletop positions. Positions are adjusted with remote control, in a quite and nonvibrating manner at constant speed.

Surgiline 2000E Surgical Table is equipped with battery recharge system. Movements of the table are fed from battery circuit. Battery power of the table has a 30 minutes operating time in cases of mains failure providing approximately 30 complete position changes.

The tabletop consists of four sections: head, backrest, seat and leg, each with removable, easy to clean, anti-bacterial, and waterproof 80 mm mattress for patient comfort.

Easy to install radiographic tabletop kit\* enables using Surgiline 2000E for X-ray procedures. Sliding motion of the tabletop enables Surgiline 2000E Surgical Table compatible for use with C-Arm radiological equipment.

The durable cast iron table base and electrohydraulic piston is covered with stainless steel. The "T" shape design of the Surgiline 2000E allows easy access to the tabletop in any position. Special brakes and anti-static casters allow ease of mobility. Stainless steel side rails accept a variety of available accessories.

Surgiline 2000E Surgical Table is a state-of-the-art product and with both quality and competitive pricing it is an excellent choice for today's demanding surgery theaters.

(\* Optional)

### Technical Details

- Height adjustment:** 725 - 1055 mm
- Height adjustment (w/mattress):** 805 - 1135 mm
- Tabletop:** 500 x 1880 mm
- Table base:** 370 x 830 mm
- Head section:** 500 x 240 mm
- Backrest section:** 500 x 490 mm
- Seat section:** 500 x 510 mm
- Right leg section:** 235 x 555 mm
- Left leg section:** 235 x 555 mm
- Voltage:** 230V, 50Hz
- Remote control voltage:** DC 5V
- Motor inputs:** DC 24V, 10A
- Battery:** 2 each, DC 12V, 7Ah, maintenance free, lead acid battery
- Battery recharge voltage:** DC 27.2V
- Battery recharge time:** 18h (full charge)
- Low Battery Level Alarm:** Sound alarm

### Packing Details

- Net weight:** 250 kg
- Gross weight:** 295 kg
- Packing dimensions:** 1000 x 1680 x 1100 mm

### Positions

- Head section:** Upward 15°, Downward 90°, Extendable 5 cm, Removable
- Backrest section:** Upward 75°, Downward 45°, Upward 45° (w/radiographic tabletop and mattress)
- Leg section:** Upward 15°, Downward 90°, Spreading 90°, Removable, Independent movement
- Lateral:** Right 20°, Left 20°
- Trendelenburg:** 25°
- Reverse Trendelenburg:** 25°
- Lateral nephrectomy:** 135°
- Reflex:** 105°
- Flex:** 225°
- Tabletop:** Rotates 355°
- Horizontal tabletop sliding:** 250 mm

