

Vaculine Maxi and Vaculine Maxi E Surgical Aspirator

Ref: M4E 500 0 0 - Vaculine Maxi

Ref: M4E 650 0 0 - Vaculine Maxi E

For operating room, critical care, emergency care and outpatient clinic use

- ▶ 2, 3 and 5 L single/double polycarbonate bottle combinations
- ▶ Suitable for use with canister systems
- ▶ 0-680/710 mmHg vacuum range (at sea level)
- ▶ Vacuum diversion knob to select bottle
- ▶ Vacuum adjust knob and 100mm gauge
- ▶ Output: 50L/ml Air filter
- ▶ One suction bag remains functional while other is being discarded *
- ▶ Ergonomic control panel design with 25 degree angle
- ▶ Fluid aspiration into the pump:
 - Floating valve,
 - Hydrophobic filter:
 - Tight seal lid with clear view of the filter, enabling prompt filter change
 - Hydrophobic filter allows free air-flow, but prevents fluid overflow, protecting the pump
- ▶ Steering handle
- ▶ Basket in the rear for spare suction catheters*
- ▶ Rinse solution receptacle and two separate suction line holders on the sides
- ▶ Polycarbonate bottle/canister compatible with steam sterilization (Applicable for 2 and 3 L bottles/canisters)
- ▶ Operation status indication lamp
- ▶ Silent, dry vacuum pump
- ▶ Design features to prevent airborne contamination in operating rooms and clean rooms. No requirement for periodic lubrication and maintenance
- ▶ Motor speed: 1400 rpm
- ▶ Sturdy body construct with durable polyester dye
- ▶ Mobile
- ▶ Four 360° swivel castors with 100 mm diameter, two with brakes
- ▶ Rated voltage: 230 V, 50/60 Hz
- ▶ Power consumption: 145W
- ▶ Dimension (Vaculine Maxi): 485 x 650 x 985 mm
- ▶ Dimension (Vaculine Maxi E): 485 x 655 x 965 mm
- ▶ Weight (Vaculine Maxi): 35 kg
- ▶ Weight (Vaculine Maxi E): 32 kg
- ▶ Packing dimensions: 630 x 630 x 1080 mm

Vaculine Maxi and Vaculine Maxi E Surgical Aspirator is ordered as the basic unit. The bottle kits configuration can be selected according to customer requirements and shall be ordered separately.



Vaculine Maxi



Vaculine Maxi

Vaculine Maxi E

* Available for Vaculine Maxi; not available for Vaculine Maxi E



Vaculine Maxi E

Vaculine bottle sets can be ordered as:
 2 x 2 liter, 2 x 3 liter or 2 x 5 liter separately
 according to customer requirements.



Vaculine Bottle Sets
 Ref: M4X 710 0 0 - 5 L
 Ref: M4X 703 0 0 - 3 L
 Ref: M4X 702 0 0 - 2 L

Vaculine Mini Aspirator

Portable

Ref: M8F 200 0 0

For bedside application

- ▶ 1 Liter collecting bottle
- ▶ Suitable for use with canister systems
- ▶ 0-600/620 mmHg vacuum range (at see level)
- ▶ Output: 20 L/min.
- ▶ Air filter
- ▶ Fluid aspiration into the pump:
 - Floating valve,
 - Hydrophobic filter:
 - Tight seal lid with clear view of the filter, enabling prompt filter change
 - Filter allows free air-flow, but prevents fluid overflow, protecting the pump
- ▶ A separate suction line holder on the side
- ▶ Polycarbonate canister compatible with steam sterilization
- ▶ Operation status indication lamp
- ▶ Silent, dry vacuum pump
- ▶ Design features to prevent airborne contamination in operating rooms and clean rooms. No requirement for periodic lubrication and maintenance
- ▶ Motor speed: 1400 rpm
- ▶ Mobile
- ▶ Sturdy body construct with durable polyester dye
- ▶ Rated voltage: 230V, 50/60Hz
- ▶ Power consumption: 60W
- ▶ Dimensions:
 - Device : 345 x 250 x 305 mm
 - Stand : Ø620 mm, Height: 720 mm
- ▶ Weight :
 - Device : 6 kg
 - With stand: 12.5 kg



Available with stand
Stand Ref: M8F 230 5 0



Vaculine Aspirator

Thoracic Drainage

Ref: M4Z 200 0 0

For thoracic drainage procedures

- ▶ 3 L collection bottle made of polycarbonate material
- ▶ 25 cmH₂O vacuum power (at sea level)
- ▶ Intermittent vacuum
- ▶ Overflow safety device for pump
- ▶ Water manometer with cm scaled
- ▶ There is a safety jar that prevents fluid access to the vacuum pump from the water manometer
- ▶ Fan for cooling
- ▶ Phosphating process is applied before dyeing and powder paint (epoxy polyester) is housing
- ▶ Mobile
- ▶ Supply voltage: 230V, 50/60 Hz
- ▶ Energy consumption: 350 W
- ▶ Dimensions: 485 x 660 x 965 mm
- ▶ Weight: 28 kg
- ▶ Packing dimensions: 645 x 635 x 1030 mm



Vaculine Aspirator

Lipectomy

Ref: M4T 000 0 0

For lipectomy applications

- ▶ 2 L and 1 L collecting bottles with parallel connection
- ▶ Polycarbonate canisters compatible with steam sterilization
- ▶ 0-680/710 mmHg vacuum range
- ▶ Single or double pump use option
- ▶ Pump displacement: 90 L/min with single pump, 130 L/min with double pump
- ▶ Vacuum adjustment
- ▶ Vacuum gauge
- ▶ Quiet operation
- ▶ Maintenance and oil free piston pump
- ▶ Power consumption: 480W
- ▶ Mobile
- ▶ Rated voltage: 230V, 50/60Hz
- ▶ Dimensions: 570 x 550 x 860 mm
- ▶ Weight: 78 kg
- ▶ Packing dimension: 620 x 620 x 920 mm



Vaculine Collecting Bottles and Bottle Cap



Vaculine Collecting Bottle 5 L
Ref: M4X 753 1 0



Vaculine Collecting Bottle 3 L
Ref: M4X 731 1 0



Vaculine Collecting Bottle 2 L
Ref: M4X 721 1 0



Vaculine Collecting Bottle 1 L
Ref: M4X 711 1 0



Vaculine Collecting Bottle Cap w/Hydrophobic Filter
Ref: M4X 372 5 0